

TECHNICAL NOTE

Assay Scheme and Configuration of Chromium™ Single Cell V(D)J Libraries

INTRODUCTION

The Chromium™ Single Cell V(D)J Solution produces Illumina® sequencer-ready libraries for profiling T-cell receptor (TCR) or B cell immunoglobulin (Ig) repertoires from single cells. In addition, users have the option of profiling 5' gene expression from the same cells if cDNA amplification is performed prior to target enrichment for TCR and/or Ig transcripts. During library preparation, sequence components essential for Illumina sequencing and downstream data analysis are incorporated into the final library construct. The sequence components are introduced via Gel Beads and during the library preparation steps of the workflow.

Each Gel Bead contains millions of oligo primers that are comprised of the following sequences (Figure 1):




Figure 1. Schematic of a Single Cell 5' Gel Bead oligo primer.


10x Genomics® GemCode™ Technology partitions thousands of cells into droplets, each called a Gel Bead-in-Emulsion (GEM). Once partitioned, the Gel Bead dissolves and its oligo primers are released into the aqueous environment of the GEM. The cell captured in the GEM is also lysed. The contents of the GEM (oligos, lysed cell components and Master Mix that contains the Poly-dT RT primer) are incubated in a reverse transcription (RT) reaction to generate full-length, barcoded cDNA from poly-adenylated mRNA. The reverse transcriptase incorporates the Gel Bead oligo via a template switching reaction at the 5' end of the transcript. All cDNA generated within an individual GEM share a common 10x Barcode.

The Single Cell V(D)J Solution offers the option to generate:

- i. Protocol Option 1: Direct Target Enrichment - enriched library from either T cells or B cells, directly from first-strand cDNA
- ii. Protocol Option 2: cDNA Amplification followed by Target Enrichment - enriched T cell library and/or an enriched B cell library, and/or a 5' gene expression library from amplified cDNA from the same cells

The Single Cell V(D)J Protocols produce V(D)J enriched and 5' gene expression Illumina® sequencer-ready libraries. A library comprises standard Illumina paired-end constructs which begin and end with P5 and P7, respectively. For V(D)J enriched libraries, Read 1 encodes the 16 bp 10x™ Barcode, 10 bp UMI, and 13 bp Switch Oligo, as well as the 5' end of an enriched transcript. For 5' gene expression libraries, Read 1 encodes the 16 bp 10x Barcode and 10 bp UMI. Due to Enzymatic Fragmentation, for both libraries Read 2 encodes a random internal fragment of the corresponding insert. Sample index sequences are incorporated as the i7 index read. A schematic of the final library constructs is shown in Figure 2.

V(D)J Enriched Library Structure:



5' Gene Expression Library Structure:




Figure 2. Schematic of final library constructs and recommended sequencing run parameters for the Single Cell V(D)J Protocol options.

An overview of the Single Cell V(D)J Protocol options and how individual sequence components are incorporated during library construction is presented in Figure 3 and Figure 4. Table 1 provides detailed reaction products and oligo sequences for Protocol Option 1: Direct Target Enrichment and Table 2 provides detailed reaction products and oligo sequences for Protocol Option 2: Target Enrichment from Amplified cDNA and 5' Gene Expression. Refer to the *Chromium™ Single Cell V(D)J Reagent Kits User Guide* (CG000086) for more details.




Figure 3. Assay schematic for Single Cell V(D)J Protocol Option 1: Direct Target Enrichment Library Construction. Fragmentation occurs through the length of the V(D)JC gene segments. Only fragmented products that contain the C region are shown for simplicity. TSO = Template Switch Oligo.




Figure 4. Assay schematic for Single Cell V(D)J Protocol Option 2: Target Enrichment from Amplified cDNA and 5' Gene Expression Library Construction. Fragmentation occurs through the length of the V(D)JC gene segments. Only fragmented products that contain the C region are shown for simplicity. TSO = Template Switch Oligo.

CONCLUSION

We have presented a detailed description of the assay configuration for the Single Cell V(D)J Solution, including enriched libraries and 5' gene expression libraries. Individual steps during library construction outlined here provide additional insight and may serve as a reference to customize the library preparation workflow for unsupported technical development and applications.

REFERENCES

- *Chromium™ Single Cell V(D)J Reagent Kits User Guide* (CG000086)

Protocol Step 1.5 – GEM-RT Incubation

Gel Bead Oligo Primer (TSO) (PN-220112)	
Poly-dT RT Primer (PN-2000007)	
Reverse Transcript Product	

Protocol Step 3.1 – Direct Target Enrichment 1

Human T Cell Mix 1 (PN-2000008)	<p>Forward primer: (final concentration of 2 μM) PCR Primer 5'-AATGATAACGGCGACCACCGA-GATCTACACTCTTCCTACACGACGCTC-3'</p>	<p>Reverse primers: (final concentration of 1 μM each) 5'-TGAAGGCCTTGACATGCA-3' 5'-TCAGGCAGTATCTGGAGTCATGAG-3'</p>
Human B Cell Mix 1 (PN-2000035)	<p>Forward primer: (final concentration of 1 μM) PCR Primer 5'-AATGATAACGGCGACCACCGA-GATCTACACTCTTCCTACACGACGCTC-3'</p>	<p>Reverse primers: (final concentration of 0.5 μM each) 5'-CAGGGCACAGTCACATCCT-3' 5'-TGCTGGACACGCATTGTA-3' 5'-GGTTTGTGTCGACCCAGTCT-3' 5'-TTGTCCACCTTGGTGTGCT-3' 5'-CATGACGCTTGGAAAGGCA-3' 5'-TGTGGACTTCCACTG-3' 5'-TTCTCGTACTCTGCTTTGCTCAG-3'</p>

Protocol Step 3.3 – Direct Target Enrichment 2

Human T Cell Mix 2 (PN-2000009)	<p>Forward primer: (final concentration of 2 μM) PCR Primer 5'-AATGATAACGGCGACCACCGA-GATCT-3'</p>	<p>Reverse primers: (final concentration of 1 μM each) 5'-AGTCCTCAGCTGGTACACG-3' 5'-TCTGATGGCTCAAACACAGC-3'</p>
Human B Cell Mix 2 (PN-2000036)	<p>Forward primer: (final concentration of 1 μM) PCR Primer 5'-AATGATAACGGCGACCACCGA-GATCT-3'</p>	<p>Reverse primers: (final concentration of 0.5 μM each) 5'-GGGAAGTTCTGGCGGTCA-3' 5'-GGTGGTACCCAGTTATCAACGAT-3' 5'-GTGTCACAGTCACCATCAC-3' 5'-TCCTGAGGACTGTAGGACAGC-3' 5'-CACGCTGCTCGTATCCGA-3' 5'-TAGCTGCTGGCCGC-3' 5'-GCGTTATCCACCTTCACTGT-3'</p>

Targeted Amplification Product	
--------------------------------	--

Protocol 7.2 – Adaptor Ligation

Table 1. Detailed reaction products and oligo sequences for the Single Cell V(D)J Protocol Option 1: Direct Target Enrichment Library Construction. Protocol steps correspond to the *Chromium™ Single Cell V(D)J Reagent Kits User Guide* (CG000086).

Protocol Step 1.5 – GEM-RT Incubation

Protocol 4.1 – cDNA Amplification

Protocol Step 5.1 – Target Enrichment 1 from Amplified cDNA

Human T Cell Mix 1 (PN-2000008)	Forward primer: (final concentration of 2 μ M) PCR Primer 5' -AATGATACGGCACCACCGA-GATCTACACTTTCCACACGACGCTC-3' 	Reverse primers: (final concentration of 1 μ M each) 5' -TGAAGCGTTGCACATGCA-3' 5' -TCAGGCACTGAGTCATTGAG-3' 	Enrichment Outer Primer
Human B Cell Mix 1 (PN-2000035)	Forward primer: (final concentration of 1 μ M) PCR Primer 5' -AATGATACGGCACCACCGA-GATCTACACTTTCCACACGACGCTC-3' 	Reverse primers: (final concentration of 0.5 μ M each) 5' -CAGGGCACAGTCACATCCT-3' 5' -TGTGGACCACGCCATTGTA-3' 5' -GGTTTTGTCGACCCAGTCT-3' 5' -TTGCCACCTTGGTGTGCT-3' 5' -CATGACGTCTTGAAGGCA-3' 5' -TGIGGGACTTCACTG-3' 5' -TTCTCGTAGTCTGTTGCTCAG-3' 	Enrichment Outer Primer

Protocol Step 5.3 – Target Enrichment 2 from Amplified cDNA

Human T Cell Mix 2 (PN-2000009)	Forward primer: (final concentration of 2 μ M) PCR Primer 5' - AATGATAACGGCACCACCGA-GATCT-3'	Reverse primers: (final concentration of 1 μ M each) 5' - AGTCTCTCAGCTGGTACACG-3' 5' - TCTGATGGCTCAAACACAGC-3'	Enrichment Inner Primer
Human B Cell Mix 2 (PN-2000036)	Forward primer: (final concentration of 1 μ M) PCR Primer 5' - AATGATAACGGCACCACCGA-GATCT-3'	Reverse primers: (final concentration of 0.5 μ M each) 5' - GGGAAGTTCTGGCGGTCA-3' 5' - GGTGGTACCCAGTTATCAAGCAT-3' 5' - GTGTCCCAGGTACACATCAC-3' 5' - TCCTGAGGACTGTAGGACAGC-3' 5' - CACGCTGCTCGTATCCGA-3' 5' - TAGCTGCTGGCCGC-3' 5' - GCGTTATCCACCTTCACTGT-3'	Enrichment Inner Primer

Targeted Amplification Product	<p>5' - AATGATAACGGCACCACCGA - GATCTACACTCTTCCCTACACGACGCTCTCCGATCT - NNNNNNNNNNNNNNNN - NNNNNNNNNNNN - TTTCTTATATGGG - cDNA_Insert - Inner_Primer-3' 3' - TTACTATGCCGCTGGTGGCT - CTAGATGTGAGAAAGGGATGTGCTCGAGAACGGCTAGA - NNNNNNNNNNNNNNNN - NNNNNNNNNNNN - AAAGAATATACCC - cDNA_Insert - Inner_Primer-5'</p>
	Protocol 7.2 – Adaptor Ligation
Adaptor (Read 2) (PN-220026)	<p>5' - GATCGGAAGAGCACACGCTGAACTCCAGTCAC - 3' 3' - TCTAGCCTTCTCG - 5'</p>
	Ligation Product
Sample Index PCR Primer (PN-220111 and PN-220103)	<p>Forward primer: SI-PCR Primer P5 – Partial Read 1 5' - AATGATAACGGCACCACCGA - GATCTACACTCTTCCCTACACGACGCTC - 3'</p> <p>Reverse primer: Chromium™ i7 Sample Index P7 – Sample Index – Partial Read 2 5' - CAAGCAGAAGACGGCATACGAGAT - NNNNNNNN - GTGACTGGAGTTCAAGCTGT - 3'</p>
	Sample Index PCR Product

Table 2A. Detailed reaction products and oligo sequences for the Single Cell V(D)J Protocol Option 2: Target Enrichment from Amplified cDNA Library Construction. Protocol steps correspond to the *Chromium™ Single Cell V(D)J Reagent Kits User Guide* (CG000086).

Protocol Step 1.5 – GEM-RT Incubation

Gel Bead Oligo Primer (TSO) (PN-220112)	
Poly-dT RT Primer (PN-200007)	
Reverse Transcript Product	
Protocol 4.1 – cDNA Amplification	
cDNA Primer Mix (PN-220106)	
cDNA Amplified Product	
Protocol 6.3 – Adaptor Ligation	
Adaptor (Read 2) (PN-220026)	
Ligation Product	
Protocol 6.5 – Sample Index PCR	
Sample Index PCR Primer (PN-220111 and PN-220103)	
Sample Index PCR Product	

Table 2B. Detailed reaction products and oligo sequences for the Single Cell V(D)J Protocol Option 2: 5' Gene Expression Library Construction. Protocol steps correspond to the *Chromium™ Single Cell V(D)J Reagent Kits User Guide* (CG000086).

Notices

Document Number

CG000109 Rev D *Technical Note*

Legal Notices

© 2017 10x Genomics, Inc. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, Inc., is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. "10x", "10x Genomics", "Changing the Definition of Sequencing", "Chromium", "GemCode", "Loupe", "Long Ranger", "Cell Ranger" and "Supernova" are trademarks of 10x Genomics, Inc. All other trademarks are the property of their respective owners. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Product(s) in practicing the methods set forth herein has not been validated by 10x, and such non-validated use is NOT COVERED BY 10X STANDARD WARRANTY, AND 10X HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE.

Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics, Inc., terms and conditions of sale for the Chromium™ Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics, Inc that it currently or will at any time in the future offer or in any way support any application set forth herein.

Customer Information and Feedback

For technical information or advice, please contact our Customer Technical Support Division online at any time.

Email: support@10xgenomics.com

10x Genomics 7068 Koll Center Parkway
Suite 401
Pleasanton, CA 94566 USA